

Image Analysis Basics

Introduction to the Digital Image

Pol Biolimage Analysis Symposium

- Training School -

Jan Brocher

Acknowledgement

Robert Haase

Marcelo Leomil Zoccoler

PoL for hosting

You for participating, your interest and questions !



Let's get started...

...at the beginning



What is a Pixel




Is a Pixel a Square ?







A Pixel IS NOT a SQUARE!

Pixel = picture element = point-sample of a real world signal




Alternative Pixel Representations – There is no Optimum







From Points to Pixels




Fluorescent
Molecules
(e.g. beads)



Occurrence
in imaging system
as PSF




Areas covered by
digital pixel representations




digital pixels

Pixels = Data !



Always save your Raw Data in the Original File Format



Original
(commercial)
File type

Metadata = **368 lines**




Exported
Tiff

Metadata = **28 lines**


92% metadata loss

DON'T USE JPEG

Original




JPEG (66% compression)



Metadata = 0 lines


Pixel Gray Value reflect the Recorded Intensity



7	6	19	45	85	99	95	85	84	75	75	78
4	5	14	45	79	83	87	94	76	66	80	
4	6	13	40	86	91	84	94	104	77	64	
4	6	9	42	88	119	139	110	105	93	75	
4	4	6	26	99	166	236	243	148	101	96	
4	4	3	15	87	194	255	255	230	100	93	
3	3	4	7	44	117	238	253	178	102	97	
4	3	4	5	36	60	111	132	100	74	75	
3	5	5	6	25	46	65	73	74	71	62	
6	3	4	5	25	45	74	74	82	71	79	
6	5	4	5	17	54	65	94	88	87	88	

Counting all Pixels of each Intensity = Histogram


7	6	19	45	85	99	95	85	84	75	78
4	5	14	45	79	83	87	94	76	66	80
4	6	13	40	86	91	84	94	104	77	64
4	6	9	42	88	119	139	110	105	93	75
4	4	6	26	99	166	236	243	148	101	96
4	4	3	15	87	194	255	255	230	100	93
3	3	4	7	44	117	238	253	178	102	97
4	3	4	5	36	60	111	132	100	74	75
3	5	5	6	25	46	65	73	74	71	62
6	3	4	5	25	45	74	74	82	71	79
6	5	4	5	17	54	65	94	88	87	88





Dynamic / Intensity range

Counting all Pixels of each Intensity = Histogram


7	6	19	45	85	99	95	85	84	75	78
4	5	14	45	79	83	87	94	76	66	80
4	6	13	40	86	91	84	94	104	77	64
4	6	9	42	88	119	139	110	105	93	75
4	4	6	26	99	166	236	243	148	101	96
4	4	3	15	87	194	255	255	230	100	93
3	3	4	7	44	117	238	253	178	102	97
4	3	4	5	36	60	111	132	100	74	75
3	5	5	6	25	46	65	73	74	71	62
6	3	4	5	25	45	74	74	82	71	79
6	5	4	5	17	54	65	94	88	87	88



Histogram Inspection reveals Usage of the Dynamic Range






Histogram Inspection reveals Changes in Images





applied
contrast
enhancement

(point operation)



Bit-Depth defines Range of Image Intensity Values


- A bit is the smallest memory unit in computers, *atomic data*.
- The bit-depth n enumerates how many different intensity values are present in an image:
 - 2^n grey values
- In microscopy, images are usually stored as 8, 12, or 16-bit images.







Changes in Bit-Depth change the Pixel Intensity Range



\neq



\neq







Changes in Bit-Depth change the Pixel Intensity Range



≠








≠



Examples of Intensity-based versus True Color Images



Intensity-based images

8-, 10-, 12-, 14-, 16- and 32-bit





True-color images

24- / 48-bit






True-color channels

8-bit



Colormaps / Pseudocolor / Lookup tables (LUTs)



- The lookup table decides how the image is displayed on screen.
- Applying a different lookup table does not change the pixel values, just there appearance



Colormaps / lookup tables



- Choose visualization of your color tables wisely!
- Think of people with red/green blindness!

Common view









Replace red
with
magenta!




Red/green blind people
may see it like this





Choose LUTs / Colormaps wisely !



Usability of LUTs / Colormaps depends on Context



LUTs can be used to detect Uneven Lighting







Uneven illumination



Uneven illumination
needs to be corrected
before analysis

Why do we need to address illumination problems?

4 objects of the same size, **NO** uneven illumination



4 objects of the same size, **WITH** uneven illumination




There is no such thing as “Background-Noise”

“Background-Noise” ≈




Sometimes you don't see your signal in all the other signals




Unspecific signals need correction to not influence image analysis negatively

Difference between Background and Noise and Unspecific staining








You can Subtract Background but not Noise






Background
(low noise)




Strong Noise
(and background)



Noise Reduces the Precision of your Analysis




Precision and Accuracy in Quantitative Measurements



Effect of Background and Noise on Analysis Precision and Accuracy

high noise



high background

